

How to use the DryKids Wireless bed wetting alarm

- This new alarm is a game-changer, making it possible to wake children in several ways.
- If the child is under age 6, the receiver can be in the parents room.
- If the child had taken an alarm off previously, the receiver can be placed across the bedroom.
- If the child tends to remove a clip, the sensor can be covered with another pair of briefs
- If there is a problem with placement, it can be clipped to a mini-pad in the underwear

See below for more tips...

Instructions for use:

First, plug the charging cable into one end of the case. When both LED lamps are green, it is charged.

Remove the clear film off the receiver's screen.

Turn the receiver on with a switch on the back.

Choose the Sound setting by pressing the 'set' button several times, then hit "ok". You can also choose vibration, or both sound and vibration. NOTE that vibration is rarely used, because it won't wake the child nor cause the 'aversion' reaction. After a week or two the child will hear the alarm or wake before it sounds and get up. Vibration is too easy to ignore.

Clip the sensor to the side or fly of 'tighty whities' where it will get wet. Note that the contacts are on the inside of clip, on the larger side. If applying the sensor to a mini-pad, the larger part of sensor should be on the inside, where wetness occurs.

Silence the alarm by turning it off with the switchcn.



We found recently that it's easy to forget to turn on or off the alarm, and it does like to be fully charged. So best idea is to keep it in the case at bedside, always on and always charging.